
SUDORmed 1200

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User manual

Mode d'emploi

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
Deutsch.....	2
English.....	14
Français.....	26
Nederlands.....	39

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Iontophorese-System

SUDORmed 1200

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 **Leitungswasser-Iontophorese**

 **Medikamenten-Iontophorese**

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
 *Tap water iontophoresis*

Tables of Content

Summary.....	15
Before operation	15
Safety instructions - side effects - contraindications	15
Further information for your safety.....	16
Intended usage, effectiveness and treatment plan	17
Tap water iontophoresis for the treatment of hyperhidrosis („excessiv sudation“).....	17
Drug iontophoresis for infiltration of drugs through the skin	17
Components of the Iontophoresis-Systems SUDORmed 1200	18
Control unit of the SUDORmed 1200 with LCD Display	18
Accessories	19
Starting up the Iontophoresis-System SUDORmed 1200 and therapy performance	19
Starting up and therapy performance of the hyperhidrosis treatment	19
Starting up and therapy performance of the drug iontophoresis.	20
Application of the therapy accumulator.	22
Rechargement or replacement of the accumulators.....	22
Care, maintenance and disposal.....	23
Security Technical Control (STK).....	23
Guarantee	23
Technical specifications.....	23
Sign and symbols	24
Frequently asked questions	25

Before operation

Dear customer,

We are happy that you bought this product from DAVITA® and congratulate you on your choice.

You have bought a product which can be a contribution for the improvement of your quality of life for many years. In this user manual particularly the functions of the iontophoresis system are described in detail.

Before using the device, the user manual should be read carefully.

Safety instructions – side effects – contraindications

Up to now there are noted no dangerous side effects.

Attention:

It is not allowed to use the iontophoresis system regarding patients:

- ① with a heart pacemaker
- ① with an ICD (implanted cardiofibrillator)
- ① during pregnancy

- ① with metalliferous Intrauterine devices (Coil)
- ① with metal implants in the area of current flow (arms or legs)
- ① with large skin defects / wounds, which can not be masked with Vaseline and existing scar tissue
- ① in the temporal regions and orbital region
- ① with tumour diseases
- ① with thrombosis
- ① with allergies or sensitivities to the drugs, which are applied
- ① with strongly reduced sensibility of hands and feet (e.g. polyneuropathy)
- ① Children younger than 5 years

Warning:

- ① Iontophoresis treatment can cause skin irritations and burns. It is important to note the following:
 - 1) Direct current used in iontophoresis can cause erythema which will generally resolve within a few hours to a few days.
 - 2) The user should report any undue burning or pain during treatment at once. Pause the treatment, inspect the area under the electrodes and make any necessary corrective actions before resuming the treatment.
 - 3) Before the treatment, remove metallic jewellery (rings, bracelets etc). If metallic jewellery comes into contact with the water bath the electrical current concentration at these points can lead to light burns.
- ① Failure to observe the following precautions may result in excessive skin irritations or burns:
 - 1) Electrodes should not be pressed during the treatment
 - 2) Do not use electrodes which have been altered or appear damaged
 - 3) Do not apply electrodes over damaged skin.
- ① Exercise caution in handling the unit. It is not allowed to drop or immerse the control unit into water. It is not allowed to connect the unit to external devices. Doing so may cause a malfunction or patient injury.
- ① It is not allowed to apply electrodes over or across the right and left temporal regions (area of the temples), or use the SUDORmed 1200 and electrodes for treatment of the orbital region (area of the eyes). Otherwise it can lead to transient visual disturbances.
- ① Users with known sensitivity to electrical current should be treated with lower current settings than those recommended for general use. If a treatment results in prolonged skin irritations or burns, additional treatments should not be done and a physician has to be consulted.
- ① Patients should be asked about their history of drug allergies or sensitivities. The SUDORmed 1200 and the electrodes should not be used on any patient who has a known allergy or sensitivity to the drug being administered. The drug package insert has to be consulted for additional contraindications and warnings.
- ① Only use the non-conductive plastic trays provided.
- ① It is not allowed to use two units at the same time while treating one patient.

Further information for your safety

- ✚ Place the therapy unit on a flat surface.
- ✚ The device should be at room temperature when it is switched on.
- ✚ The unit should not be used near by short waves- and microwave therapy devices, Because it could result in making the iontophoresis device provide an incorrect electrical current. The distance should be at a minimum 1 meter.
- ✚ A treatment voltage of more than 2 volt requires an increased attention. The electrodes should not put on selectively because it can lead to very high current density and can trigger skin irritations at these places.
- ✚ Portable and mobile high frequency equipment (e.g. mobile telephones) can affect the unit.
- ✚ You must avoid using accessories which are not approved by the manufacturer because it could lead to a damage of the unit and unsafe use. Also it can lead to an increased emission of electromagnetic waves or to a reduced interference resistance.
- ✚ Don't open the unit. All servicing must be carried out by a company which is authorised by DAVITA®
- ✚ This device must only be used indoors. Do not use in a room with high humidity. Keep water away from the Control unit.

Intended usage, effectiveness and treatment plan

The DAVITA® Iontophoresis-System SUDORmed 1200 can be optimally used for tap water iontophoresis and drug iontophoresis.

Tap water iontophoresis for treatment of hyperhidrosis („excessive sudation“)

- ✚ The DAVITA® Iontophoresis-System SUDORmed 1200 is used for the treatment of hyperhidrosis (excessive sweating) of the hands, feet and underarms.
- ✚ During the iontophoresis-therapy a low and safe therapy current is conducted through the affected parts of the body. This current flow is enabled because of the water bath and results in decreased sweating of the parts of the body which are in the therapy water.
- ✚ This therapeutic effect has been verified in numerous medical studies. It is not entirely understood how iontophoresis works. Physicians believe that the current flow irritates the synoptical transitions of the sweat gland nerves so strongly, that they are not stimulated to secrete sweat. This means that the “nervous supply line“ is affected and not the sweat gland itself.
- ✚ The strength of the treatment voltage can be controlled to your personal preference, but will not exceed maximum defined values, ensuring safe usage.

The iontophoresis-therapy can be classified into two phases:

- Phase 1: During phase 1, you have to carry out one treatment per day for 20 minutes. After circa 10 – 15 treatments the sweating will be normalized. One treatment per day can be done for a period of 2 – 3 months to reach a permanent relief.
- Phase 2: Successful Iontophoresis therapy requires periodic treatments over the years in order to keep sweating normalised. Depending on the severity of the patient's hyperhidrosis, sweating will continue to be normalised with a further one to three treatments per week for approx. 20 minutes.


Drug iontophoresis for the infiltration of drugs through the skin

- ✚ The iontophoresis system SUDORmed 1200 can also be used for drug iontophoresis ions walk from the + to the – pole or inverse. Thereby drugs are infiltrated through the skin into the tissue joint and into the bloodstream. This allows for example to increase considerably the positive effect of the anti-inflammatory and analgetic of ointments and gels. The application occurs for example by rheumatism, bruise, strains and arthritis.
- ✚ On the package insert of the drug it must be stated if the pharmaceutical product is suitable to iontophoresis. In the package insert should also be specified, on which electrode the drug must be applied. At minus (-) apply the drug to the electrode, which is connected with the white cable and the red plug. Plus (+) on the electrode, which is connected with the white cable and the white plug.
- ✚ The cables have to be inserted according to their colours into the sockets of the device. If not specified, you take the drug on both electrodes and change the polarity after half of the treatment time.
- ✚ The treatment voltage must be set very fine and should not exceed 5 volt.
- ✚ The treatment time takes about 20 minutes, daily.

Components of the Iontophoresis System SUDORmed 1200


Your DAVITA®-therapy system consists of a control unit integrated in the case and further accessories, which you will get to know in the following chapter.

Control unit of the SUDORmed 1200 with LCD Display


Functions of button

SUDORmed 1200 has the feature to control channel A and channel B individually. Therefore there are 2 multifunction buttons for the A and B channel.

- ✚ Multifunction button A controls channel A. Multifunction button B controls channel B. For saving energy, while setting the multifunction button A, only LCD display for A channel bright, B channel is invisible.
- ✚ Always use the ↑ or ↓ to locate on the different options of mem, timer and V (treatment voltage)
- ✚ Always use + or – to increase or decrease the setting of option.

LCD Display

Accessories

Iontophoresis-System SUDORmed 1200, Art. No.: 30035, consisting of:

- 1 x Control unit SUDORmed 1200, Art. No.: 30049
- 2 x large flexible electrodes, 140 x 200 mm, Art. No.: 35003
- 1 x connection cable, 2 m, Art. No.: 35005
- 4 x Accumulators (Type Ni-MH AA 1,2V)
- 1 x Charging device, Art. No.: 35006
- 1 x Storage box for control unit, connection cable, accumulators and charging device, Art. No.: 00475
- 1 x Case out of plastic with 2 separable trays, Art. No.: 00507
- 1 x multilingual user manual, Art. No.: 00473

Available as accessory are:

- 2 x axles electrodes with sponge case, Art. No.: 35009 for application under the axles
- 2 x small flexible electrodes, 55 x 85 mm each with sponge case, Art. No.: 35004 for the drug-iontophoresis or for application by hyperhidrosis under the axles

Starting up the Iontophoresis system SUDORmed 1200 and therapy performance

In this chapter the operation mode of the system is described in detail.

Starting up and therapy performance of the Hyperhidrosis treatment

A. For beginning the therapy all metallic items (e.g. jewellery) must be removed. **For the hyperhidrosis treatment** the red and white plug of the white connecting cable have to be connected with the large flexible electrodes. It is important to ensure that the plug is fully inserted into the electrode. The metal part of the pin should not be seen. **The large electrodes must be laid each into an tray with the black side on the button. Only the grey side of the electrodes is allowed to touch the hands or feet.** Both trays have to be filled with lukewarm water, so that the hands are covered with water completely.

To fill and drain the trays easily it is possible to separate the trays of the case. Therefore the case has to be opened completely and the parts of the trays have to be pulled apart contrary. For example therefore the left tray has to be fixed and the right tray has to be removed downwards. Every time the trays can be connected again if it is necessary.

B. The large white plug of the connecting cable has to be connected on the top of the control unit of the SUDORMed 1200 with the plug socket of channel A or B.

C. Power on: Push the ON/OFF button to turn the power on.

D. Selecting mode:

Pulsed current: Adjust Channel A buttons
Direct current Adjust Channel B buttons

Which type of current is the better one for you, has to be decided by you or your physician. Direct current should be more effective. Many patients feel the pulsed current as more comfortable.

E. Setting mem:

Push + / - button to choose:

“00” for not memory the current used parameter
“01” for memory the current used parameter (last mem call function)

F. Setting Timer:

Push + / - button to set timer between 0 and 99 minutes.
The recommended therapy duration is 20 minutes per day. (See the additional explanations on page 4f.)

G. Setting output intensity:

Push + / - button to set the output intensity.

- ⚠ For the hyperhidrosis therapy there are available a pulsed current between 0-60 volt and a direct current between 0-30 volt. Excessive currents do not lead to a better treatment success. On the contrary they can cause skin irritations.
- ⚠ Choose the intensity of the treatment voltage so that direct current treatments can not be felt and pulsed current treatments can be felt slightly. The treatment voltage is always set to your own feelings. For each treatment, the treatment voltage is felt differently.

H. Starting treatment:

To start the treatment both hands or feet have to be placed on the grey side of the two electrodes which are placed in the trays. Now the circuit is closed and the unit increases the intensity slowly from 0 volt to the set value.

The hands can always be removed from the trays, to readjust the current or to interrupt the therapy. Therefore the unit decreases the intensity automatically to 0 volt. Thereby it can occur perfectly safe and easily noticeable static discharges.

To treat hands and feet at the same time, both feet have to be put in one tray and both hands in the other tray.

To carry out a **treatment under the axles** special axle electrodes, Art. No.: 35009 or smaller electrodes with sponge

cases, Art. No.: 35004 must be used. They are not included in the scope of delivery. The axles electrodes have the form of a sponge and can be fixed safe and easy under the axles. For fixing the smaller electrodes with sponge cases it is the best to wear a tight T-shirt under which the electrodes can be moved. **The sponge cases have to be wetted very well. Otherwise, it is possible that the intensity can not be adjusted at the unit or the current can be felt strongly as biting or different strong. The following procedure should be observed:**

The jet of tap water should flow between the white foam cushions, so that the entire sponge gets wet and saturated with water. A lightly pressure on the sponge can distribute the moisture within the sponge and remove excess water. Make sure that the sponge remains in "good wet" condition. This process should be repeated several times during the application under the axles, in order to ensure a continuous current flow is given over the entire skin area of the armpit.

I. Finishing of treatment:

After the ending of the treatment time the current flow is interrupted automatically. Hands and feet are removed from the trays and the device is switched off by pressing the ON / OFF button.

To increase the success of therapy it is recommended to change the polarity for each treatment. The change of polarity can be realized by changing the position of the red and the white plug. Often it is also recommended for particularly hard water to put some salt into the water and stir well. Thus, the conductivity of the water is improved. The addition of 1 up to several teaspoons of salt can improve the treatment success even further.

Starting up and therapy performance of the drug iontophoresis

A. Before the therapy all the metal parts (such as jewellery) must be removed.

By the drug iontophoresis the control unit must be connected with the help of the two electrodes. Therefore the red and white plug must be connected with one electrode. It is important to ensure that the plug is fully inserted into the electrode. The metal part of the pin should not be seen. The large white plug of the connecting cable has to be connected on the top of the control unit of the SUDORmed 1200 with the plug socket of channel A or B.

The electrodes are placed over the parts of the body which must be treated.

On the package insert of the drug it must be stated that the pharmaceutical product is suitable for iontophoresis. The leaflet should also specify on which electrode the product must be applied. At minus (-) take the product to the electrode which must be connected with the white plug and plus (+) on the electrode which must be connected with the red plug. If not specified, you take the drug on both electrodes and change the polarity after half of the treatment time.

What quantity of the drug is to be used, must be inquired at the pharmacy or manufacturer of the drug.

B. Turning on the device: The ON/OFF-Button must be pressed, to turn on the device.

C. Selection of operation:

Pulsed current:: Set the buttons of channel A

Direct current:: Setting the buttons of the channel B

What type of current is the better one for you, must be decided by you or your physician. For the drug iontophoresis it can be used the pulsed or the direct current. The direct current is much more effective.

D. Setting the memory function:

With the + or - button the two memory functions can be selected:

"00" the last used parameters are not saved.

"01" the last used parameters are saved and are called again by the next application.

E. Setting the timer:

With the + / - button, the treatment time can be set between 0 - 99 minutes.

The recommended treatment duration is 20 minutes daily.

F. Set the intensity:

With the + / - button the intensity can be adjusted.

For the drug iontophoresis the treatment voltage has to be adjusted very finely and should not exceed 5 volt.

G. Start of treatment:

The treatment starts as soon as the cables on channel A or channel B are connected with the device.

H. Termination of treatment:

After the set treatment time the current flow is interrupted automatically. The device is turned off by the ON / OFF button and the electrodes are removed from the respective skin areas.

Application of the therapy accumulator

The device has a therapy accumulator. To recall the up to now accumulator therapy duration proceed as follows: In channel A or B on the multi-function buttons the cursor is to press up or down until the value of the bottom of the screen next to "Timer" is flashing. Then the + button on the multifunction button of the channel which is now active is to press first and in addition even the second + button on the second multifunction. Both + buttons have to be hold for 20 seconds. After the 20 seconds on the display appear, the cumulative treatment durations for channel A and channel B. The value 06 23 for example stands for a cumulative treatment duration of 62 hours and 30 minutes.

To see the normal setting on the display, the unit is to switch on and off again with the ON / Off button with the On / Off switch is off and on again briefly to turn on.

To reset the indicator of the cumulative duration treatment to 0 instead of the two + buttons both – buttons must be hold simultaneously for 20 seconds.

Rechargement or replacement of the accumulators

-  The charging of the accumulators may only be supplied with the original charger.
-  Once the charger is connected to the control unit, it is not possible for safety reasons to turn on the unit and perform a treatment.
-  **The charging time is approx. 12 – 14 hours.**
-  The unit can also be used with normal batteries, but it is not allowed to charge the device with normal batteries. Ensure that the device is only charged with accumulators!
-  The device uses 4 AA rechargeable accumulators.

A. An accumulator charge indicator on the display lights up when the accumulators are nearly discharged. Please only use the included charger for recharging the accumulators. The connection for the charger is located on the left side of the control unit.

Note: Never use not rechargeable batteries using charger for charging.

B. The accumulator compartment is located on the back of the unit. To open the accumulator compartment cover is pressed down slightly and withdraw to bottom. There are 4 accumulators visible. These must be charged with the supplied charger before the first treatment. The polarity symbols (+) and (-) are indicated on the inside of the compartment. When the accumulators are inserted incorrectly with regard to polarity, the device will not work and it is possible that the accumulators will be destroyed. Please make sure that the lid of the accumulator compartment is completely closed after the accumulators have been installed.

Push this convex and slide downward.

C. Charging function

Note! Only use the included Charger

Accumulator type: Ni-MH 1.2V AA x 4 rechargeable accumulators

Note: Do not use not rechargeable batteries while using charger for charging

D. Security Note:

Always monitor the charging process. Accumulators never leave unguarded, or in the unit overnight. There is a risk of leakage in case of overheating. Only use accumulators of the same type and with the same capacity. Please note the temperature! It is not allowed that the water gets warmer than lukewarm.

E. Accumulator information:

NiMH accumulators reach their full capacity only after several charge - discharge cycles. After about 10 charge cycles the accumulators should be discharged completely and then reloaded. Accumulators have a self discharging effect. Therefore the accumulators need to be reloaded after long periods of non-use. **In case of non-use the accumulators have to be charged every 2 months at least.** Otherwise it is possible that the accumulators are discharged too far, that a charge is no longer possible.

Care, maintenance and disposal

- ✚ Before cleaning the device, remove the plug from the power socket and switch off the unit. For cleaning use a soft and slightly moistened cloth.
- ✚ Don't use petroleum, thinner, alcohol, wax remover or other solvents.
- ✚ Cables must always be removed by holding the connectors. If the cables become badly bent or breaks in the insulation appear, then you must replace them. The cables should not be subjected to heat or chemicals.
- ✚ The electrodes and sponges must be cleaned in lukewarm water after every use. Do not use caustic lotions.
- ✚ If you don't use the device for a long time, please take out the rechargeable batteries and store them separately.
- ✚ A special maintenance is not necessary for this device. Only DAVITA[®] or qualified companies are allowed to do repairs.
- ✚ Packaging: DAVITA[®] joints to the system „Grüner Punkt“. All parts of the packaging are recyclable materials and environmentally acceptable to dispose. LEDs and units which are signed with the crossed wheelie bin symbol are not allowed to be thrown into the domestic waste. There is the commitment to dispose such electro- and electronic "Altgeräte" (old device) separately. The communal public authority which is responsible for the disposal gives information regarding the possibilities of a regulated disposal.
- ✚ For the return and disposal of accumulators and battery the battery law is valid. National wide there are re-sampling points. The accumulators and batteries must be thrown into the green collection boxes.

Security Technical Control (STK)

With the implementation of the security technical control (STK) are only instructed competent companies / laboratories e.g. which can carry out an audit for medical products. Appropriate measuring and testing facilities are required. It is only allowed to operate the medical device with the accessories which are listed in the instruction. According to MP BetriebV for Iontophorese-systems a STK is necessary for at least every two years and after any necessary repair or processing.

Guarantee

- ✚ DAVITA[®] guarantees for this instrument a period of 24 months after date of sale. The guarantee is valid for material defects and manufacturing faults and includes the provision of spare parts as well as labour. Guarantee expires if

the instrument has been damaged, or used improperly or unauthorized intervention was carried out. The guarantee becomes invalid if the unit has not been used in accordance with the instructions.

- ✚ Keep the purchase receipt as proof and submit it together with your claim for guarantee. DAVITA® will only fulfil the guarantee upon receipt of your purchase proof.
- ✚ Consumables: electrodes, sponge cases, electrode cables and the rechargeable accumulators are excluded from the guarantee.

Technical specifications

Control unit / charging device SUDORmed 1200

Dimensions:	130 x 72 x 30 mm	(B x H x T)
Weight::	265 g	
Input:	Supply voltage:	4,8 V
	Max. power input:	500 mA
Ambient temperature	Charge at:	+ 10°C until +30°C
Output direct current	Treatment voltage:	0 – 30 volt
	Output voltage:	Max. 30 V at 2000 Ohm
Output pulsed current	Treatment voltage:	0 – 60 volt
	Output voltage:	Max. 60 V at 2000 Ohm
AC/DC Adapter/charging device	Input:	AC 110 – 240 V, 50/60 Hz
	Output:	7,5 V

Signs and symbols

The following symbols appear in this manual, on the device or on it's accessories:

Caution: Consult accompanying documents

CE Mark: conforms to essential requirements of the Medical Device Directive 93/42/EEC.

Device classification type B

DISPOSAL: Do not dispose this product as unsorted municipal waste. Collection of such waste separately for special treatment is necessary. Waste Electrical Electronic Equipment (WEEE)

GRS Common collection system batteries. Batteries do not belong I household waste

Frequently asked questions

Question 1:

After starting the system, no current flow exists.

Answer :

The following measures can restore a current flow:

- ✚ To increase the conductivity of the water one to several teaspoons of salt have to be put into the water.
- ✚ The case / tray electrodes should not be pressed too hard on the tray / case ground.
- ✚ It is to check whether the plugs red / white are oxidized on the connecting cable. If necessary, the cable must be replaced (Art. No.: 35005)
- ✚ It is to check whether the tray electrodes are oxidized on the cable connection. If necessary, the cable must be replaced (Art. No.: 35005)

Question 2:

After starting the system, no current flow exists by the accessories.

Answer :

The following measures can restore a current flow:

- ✚ The axle applicators (Art. No.: 35 009) and the sponge electrodes (Art. No.: 35 004) have to be completely saturated with water to produce a current flow. To moisten the axle applicators water has to be filled in the open side of the sponge case between the foam trays. The sponge case has to be pressed with the hand several times. It is not enough to moisten the cover with water. We recommend to humidify the sponge cases 1 – 2 times during a treatment.
- ✚ It has to be checked whether the plug contacts red / white are oxidized on the connecting cable. If necessary, the cable must be replaced (Art. No.: 35005).

Question 3:

After switching on the device nothing is shown on the display.

Answer :

- ✚ It is necessary to check if the accumulators are inserted correctly. The labeling of the + and - poles has to be noted.
- ✚ It is necessary to check if the accumulators are charged.

Question 4:

During the application of the iontophoresis therapy the device switches off.

Answer :

- ✚ The state of charge of the accumulators have to be checked.

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Telephone: +49 (0) 1805 – 911 270 • Fax: +49 (0) 1805 – 911 275


Email: mail@davita.de • www.davita.de

Subject to technical changes. Rev. No.: 11061111, Art. No.: 00473

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DAVITA® Medizinische Produkte GmbH & Co. KG • Postfach 2004 • D-47518 Kleve